

KREMEN1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KREMEN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KRM1; KREMEN;

Catalog #: A10507

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83999

Isotype: IgG

Swiss Prot: Q96MU8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a high-affinity dickkopf homolog 1 (DKK1) transmembrane receptor that functionally cooperates with DKK1 to block wingless (WNT)/beta-catenin signaling. The encoded protein is a component of a membrane complex that modulates canonical WNT signaling through lipoprotein receptor-related protein 6 (LRP6). It contains extracellular kringle, WSC, and CUB domains. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene.

To place an order, please [Click HERE](#).