

MEGF10

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 122kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MEGF10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EMARDD;

Catalog #: A10508

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84466

Isotype: IgG

Swiss Prot: Q96KG7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the multiple epidermal growth factor-like domains protein family. The encoded protein plays a role in cell adhesion, motility and proliferation, and is a critical mediator of apoptotic cell phagocytosis as well as amyloid-beta peptide uptake in the brain. Expression of this gene may be associated with schizophrenia, and mutations in this gene are a cause of early-onset myopathy, areflexia, respiratory distress, and dysphagia (EMARDD) as well as congenital myopathy with minicores. Alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).