

MPP6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MPP6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VAM1; p55T; PALS2; VAM-1;

Catalog #: A10509

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51678

Isotype: IgG

Swiss Prot: Q9NZW5

Purity: Affinity purification

For research use only.

Background:

Members of the peripheral membrane-associated guanylate kinase (MAGUK) family function in tumor suppression and receptor clustering by forming multiprotein complexes containing distinct sets of transmembrane, cytoskeletal, and cytoplasmic signaling proteins. All MAGUKs contain a PDZ-SH3-GUK core and are divided into 4 subfamilies, DLG-like (see DLG1; MIM 601014), ZO1-like (see TJP1; MIM 601009), p55-like (see MPP1; MIM 305360), and LIN2-like (see CASK; MIM 300172), based on their size and the presence of additional domains. MPP6 is a member of the p55-like MAGUK subfamily (Tseng et al., 2001)

To place an order, please [Click HERE](#).