

NOMO1

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 134kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NOMO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PM5; Nomo;

Catalog #: A10510

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23420

Isotype: IgG

Swiss Prot: Q15155

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein originally thought to be related to the collagenase gene family. This gene is one of three highly similar genes in a region of duplication located on the p arm of chromosome 16. These three genes encode closely related proteins that may have the same function. The protein encoded by one of these genes has been identified as part of a protein complex that participates in the Nodal signaling pathway during vertebrate development. Mutations in ABCC6, which is located nearby, rather than mutations in this gene are associated with pseudoxanthoma elasticum (PXE).

To place an order, please [Click HERE](#).