

SIGLEC8

Reactivity: Human Mouse Rat

Tested applications: WB

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SIGLEC8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10514

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27181

Isotype: IgG

Swiss Prot: Q9NYZ4

Purity: Affinity purification

For research use only.