

C16orf54

Reactivity: Mouse Rat

Tested applications: WB

Calculated MW: 225 aa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human C16orf54

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10519

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 105369277/283897

Isotype: IgG

Swiss Prot: Q8C708

Purity: Affinity purification

For research use only.