

CABLES2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:52KD

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human CABLES2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10525

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:81928

Isotype:IgG

Swiss Prot:Q9BTV7

Purity:Affinity purification

For research use only.