

TRDMT1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRDMT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DMNT2; DNMT2; PUMET; RNMT1; MHSAILP;

Catalog #: A10535

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1787

Isotype: IgG

Swiss Prot: O14717

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein responsible for the methylation of aspartic acid transfer RNA, specifically at the cytosine-38 residue in the anticodon loop. This enzyme also possesses residual DNA-(cytosine-C5) methyltransferase activity. While similar in sequence and structure to DNA cytosine methyltransferases, this gene is distinct and highly conserved in its function among taxa.

To place an order, please [Click HERE](#).