

TRIM24

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 117kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRIM24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PTC6; TF1A; TIF1; RNF82; TIF1A; hTIF1; TIF1ALPHA;

Catalog #: A10546

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8805

Isotype: IgG

Swiss Prot: O15164

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene mediates transcriptional control by interaction with the activation function 2 (AF2) region of several nuclear receptors, including the estrogen, retinoic acid, and vitamin D3 receptors. The protein localizes to nuclear bodies and is thought to associate with chromatin and heterochromatin-associated factors. The protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains - a RING, a B-box type 1 and a B-box type 2 - and a coiled-coil region. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).