

GIPC1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GIPC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NIP; GIPC; IIP-1; TIP-2; SEMCAP; C19orf3; Hs.6454; GLUT1CBP; RGS19IP1; SYNECTIN; SYNECTIN;

Catalog #: A10550

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10755

Isotype: IgG

Swiss Prot: O14908

Purity: Affinity purification

For research use only.

Background:

GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008)

To place an order, please [Click HERE](#).