

PITX3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PITX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ASMD; ASOD; PTX3; CTPP4; CTRCT11;

Catalog #: A10569

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5309

Isotype: IgG

Swiss Prot: O75364

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family act as transcription factors. This protein is involved in lens formation during eye development. Mutations of this gene have been associated with anterior segment mesenchymal dysgenesis and congenital cataracts.

To place an order, please [Click HERE](#).