

STK24

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STK24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MST3; STK3; MST3B; STE20; HEL-S-95;

Catalog #: A10576

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8428

Isotype: IgG

Swiss Prot: Q9Y6E0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).