

PFDN6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PFDN6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HKE2; KE-2; PFD6; H2-KE2;

Catalog #: A10580

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10471

Isotype: IgG

Swiss Prot: O15212

Purity: Affinity purification

For research use only.

Background:

PFDN6 is a subunit of the heteromeric prefoldin complex that chaperones nascent actin (see MIM 102560) and alpha- and beta-tubulin (see MIM 602529 and MIM 191130, respectively) chains pending their transfer to the cytosolic chaperonin containing TCP1 (MIM 186980) (CCT) complex (Hansen et al., 1999)

To place an order, please [Click HERE](#).