

USP24

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 294kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human USP24

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10585

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23358

Isotype: IgG

Swiss Prot: Q9UPU5

Purity: Affinity purification

For research use only.

Background:

Modification of cellular proteins by ubiquitin is an essential regulatory mechanism controlled by the coordinated action of multiple ubiquitin-conjugating and deubiquitinating enzymes. USP24 belongs to a large family of cysteine proteases that function as deubiquitinating enzymes (Quesada et al., 2004)

To place an order, please [Click HERE](#).