

IL17C

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL17C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CX2; IL-17C;

Catalog #: A10587

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27189

Isotype: IgG

Swiss Prot: Q9P0M4

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a T cell-derived cytokine that shares the sequence similarity with IL17. This cytokine was reported to stimulate the release of tumor necrosis factor alpha and interleukin 1 beta from a monocytic cell line. The expression of this cytokine was found to be restricted to activated T cells.

To place an order, please [Click HERE](#).