

WDR11

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 137kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WDR11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DR11; HH14; SRI1; BRWD2; WDR15;

Catalog #: A10590

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55717

Isotype: IgG

Swiss Prot: Q9BZH6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is located in the chromosome 10q25-26 region, which is frequently deleted in gliomas and tumors of other tissues, and is disrupted by the t(10;19) translocation rearrangement in glioblastoma cells. The gene location suggests that it is a candidate gene for the tumor suppressor locus.

To place an order, please [Click HERE](#).