

SRPRB

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SRPRB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APMCF1;

Catalog #: A10591

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 58477

Isotype: IgG

Swiss Prot: Q9Y5M8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene has similarity to mouse protein which is a subunit of the signal recognition particle receptor (SR). This subunit is a transmembrane GTPase belonging to the GTPase superfamily. It anchors alpha subunit, a peripheral membrane GTPase, to the ER membrane. SR is required for the cotranslational targeting of both secretory and membrane proteins to the ER membrane.

To place an order, please [Click HERE](#).