

TBC1D15

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 79kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TBC1D15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAB7-GAP;

Catalog #: A10593

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64786

Isotype: IgG

Swiss Prot: Q8TC07

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Ras-like proteins in brain-GTPase activating protein superfamily that share a conserved Tre-2/Bub2/Cdc16 domain. The encoded protein interacts with Ras-like protein in brain 5A and may function as a regulator of intracellular trafficking. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).