

MARCH9

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MARCH9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RNF179; MARCH-IX;

Catalog #: A10596

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92979

Isotype: IgG

Swiss Prot: Q86YJ5

Purity: Affinity purification

For research use only.

Background:

MARCH9 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH enzymes add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH9 induces internalization of several membrane glycoproteins and directs them to the endosomal compartment (Bartee et al., 2004)

To place an order, please [Click HERE](#).