

VKORC1L1

Reactivity:Rat

Tested applications:WB

Calculated MW:176aa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human VKORC1L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10600

Antibody Type:

Polyclonal Antibody

Species:Rat

Gene ID:154807

Isotype:IgG

Swiss Prot:Q6TEK5

Purity:Affinity purification

For research use only.