

GPHA2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 14kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GPHA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

A2; GPA2; ZSIG51;

Catalog #: A10601

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 170589

Isotype: IgG

Swiss Prot: Q96T91

Purity: Affinity purification

For research use only.

Background:

GPHA2 is a cystine knot-forming polypeptide and a subunit of the dimeric glycoprotein hormone family (Hsu et al., 2002)

To place an order, please [Click HERE](#).