

RELA

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RELA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p65; NFkB3;

Catalog #: A10609

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 5970

Isotype: IgG

Swiss Prot: Q04206

Purity: Affinity purification

For research use only.

Background:

NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFkB1 or NFkB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFkB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).