

TUBG1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TUBG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TUBG; GCP-1; CDCBM4; TUBGCP1;

Catalog #: A10611

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 7283

Isotype: IgG

Swiss Prot: P23258

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the tubulin superfamily. The encoded protein localizes to the centrosome where it binds to microtubules as part of a complex referred to as the gamma-tubulin ring complex. The protein mediates microtubule nucleation and is required for microtubule formation and progression of the cell cycle. A pseudogene of this gene is found on chromosome 7.

To place an order, please [Click HERE](#).