

ARPC1A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ARPC1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARPC1A; Arc40; SOP2Hs; SOP2L ;

Catalog #: A1070

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10552

Isotype: IgG

Swiss Prot: Q92747

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of seven subunits of the human Arp2/3 protein complex. This subunit is a member of the SOP2 family of proteins and is most similar to the protein encoded by gene ARPC1B. The similarity between these two proteins suggests that they both may function as p41 subunit of the human Arp2/3 complex that has been implicated in the control of actin polymerization in cells. It is possible that the p41 subunit is involved in assembling and maintaining the structure of the Arp2/3 complex. Multiple versions of the p41 subunit may adapt the functions of the complex to different cell types or developmental stages. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2010]

To place an order, please [Click HERE](#).