

CSRP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP RIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:50

RIP 1:20 - 1:50

Calculated MW: 21kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CSRP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jk

Synonym:

CSRP1;CRP;CRP1;CSRP;CYRP;D1S181E;DKFZp686M148 ;

Catalog #: A1071

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1465

Isotype: IgG

Swiss Prot: P21291

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the cysteine-rich protein (CSRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010]

To place an order, please [Click HERE](#).