

RNASE13

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200 IF 1:10 - 1:100

Calculated MW: 18kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RNASE13

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RNASE13

To place an order, please [Click HERE](#).

Catalog #: A1073

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 440163

Isotype: IgG

Swiss Prot: Q5GAN3

Purity: Affinity purification

For research use only.