

GAS7

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GAS7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10733

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8522

Isotype: IgG

Swiss Prot: O60861

Purity: Affinity purification

For research use only.