

POLR3A

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 156kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human POLR3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADDH; HLD7; RPC1; RPC155; hRPC155

Catalog #: A10737

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11128

Isotype: IgG

Swiss Prot: O14802

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is the catalytic component of RNA polymerase III, which synthesizes small RNAs. The encoded protein also acts as a sensor to detect foreign DNA and trigger an innate immune response.

To place an order, please [Click HERE](#).