

## DEFB132

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 11kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human DEFB132

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

DEFB32; UNQ827

**Catalog #:** A1074

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 400830

**Isotype:** IgG

**Swiss Prot:** Q7Z7B7

**Purity:** Affinity purification

For research use only.

**Background:**

-defensins (also designated BDs, or hBDs in human) are small cationic peptides with broad-spectrum antimicrobial activity against a variety of enveloped viruses, fungi and bacteria. Produced in mucosal epithelia and neutrophils of several species, -defensins are developmentally regulated. The family of -defensin proteins share a common defensin-motif that is characterized by multiple cysteine residues and a highly conserved tertiary structure. Besides playing a significant role in host immune defense, many -defensins also are involved in sperm maturation and capacitation. -defensin 32, also known as -defensin 132, is a 95 amino acid secreted protein that most likely contains a signal peptide sequence that requires cleavage by proteolytic enzymes in order to become biologically active.

*To place an order, please [Click HERE](#).*