

IL7R

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC

Calculated MW: 459

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human IL7R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10740

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3575

Isotype: IgG

Swiss Prot: P16871

Purity: Affinity purification

For research use only.