

TIE1

Reactivity:Human

Tested applications:WB IHC ICC

Calculated MW:1138

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human TIE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10743

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7075

Isotype:IgG

Swiss Prot:P35590

Purity:Affinity purification

For research use only.