

HSP90B1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC

Calculated MW: 803

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human HSP90B1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10746

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7184

Isotype: IgG

Swiss Prot: P14625

Purity: Affinity purification

For research use only.