

DiMethyl-Histone H3-K4

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC FC

Calculated MW: 136

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human DiMethyl-Histone H3-K4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10748

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8350

Isotype: IgG

Swiss Prot: P68431

Purity: Affinity purification

For research use only.