

CALD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 70 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CALD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Catalog #: A10768

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 800

Isotype: IgG

Swiss Prot: Q05682

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).