

CALM1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50

Calculated MW: 17 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CALM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

caM; CAMI; PHKD; CPVT4; DD132; LQT14; CALML2;

Catalog #: A10769

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 801

Isotype: IgG

Swiss Prot: P62158

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the EF-hand calcium-binding protein family. It is one of three genes which encode an identical calcium binding protein which is one of the four subunits of phosphorylase kinase. Two pseudogenes have been identified on chromosome 7 and X. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).