

CD34

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50

Calculated MW: 120 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD34

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10796

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 947

Isotype: IgG

Swiss Prot: P28906

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).