

SCARB1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:200
FC 1:20 - 1:50

Calculated MW: 80 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SCARB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLA1; SRB1; CLA-1; SR-BI; CD36L1; HDLQTL6;

Catalog #: A10799

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 949

Isotype: IgG

Swiss Prot: Q8WTV0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a plasma membrane receptor for high density lipoprotein cholesterol (HDL). The encoded protein mediates cholesterol transfer to and from HDL. In addition, this protein is a receptor for hepatitis C virus glycoprotein E2. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).