

EIF4EBP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Observed MW: Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A10852

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 1978

Isotype: IgG

Swiss Prot: Q13541

Purity: Affinity purification

To place an order, please [Click HERE](#).

For research use only.