

DHCR7

Reactivity: Human Mouse

Tested applications: WB IHC IF

Calculated MW: 475

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human DHCR7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10860

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1717

Isotype: IgG

Swiss Prot: Q9UBM7

Purity: Affinity purification

For research use only.