

H3F3A

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 15 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human H3F3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

H3F3; H3.3A;

Catalog #: A10880

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 3020

Isotype: IgG

Swiss Prot: P84243

Purity: Affinity purification

For research use only.

Background:

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene contains introns and its mRNA is polyadenylated, unlike most histone genes. The protein encoded is a replication-independent member of the histone H3 family.

To place an order, please [Click HERE](#).