

ICAM1

Reactivity: Human

Tested applications: WB IHC ICC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 89 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ICAM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BB2; CD54; P3.58;

Catalog #: A10907

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 3383

Isotype: IgG

Swiss Prot: P05362

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

To place an order, please [Click HERE](#).