

IFNAR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 90/130 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IFNAR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AVP; IFRC; IFNAR; IFNBR; IFN-alpha-REC;

Catalog #: A10913

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 3454

Isotype: IgG

Swiss Prot: P17181

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.

To place an order, please [Click HERE](#).