

C19orf10

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human C19orf10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C19orf10;EUROIMAGE1875335;IL25;IL27;IL27w;R33729_1;SF20 ;

Catalog #: A1093

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56005

Isotype: IgG

Swiss Prot: Q969H8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).