

KRT17

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF FC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:200
FC 1:20 - 1:50

Calculated MW: 48 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KRT17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PC; K17; PC2; 39.1; CK-17; PCHC1;

Catalog #: A10942

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 3872

Isotype: IgG

Swiss Prot: Q04695

Purity: Affinity purification

For research use only.

Background:

This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex.

To place an order, please [Click HERE](#).