

ORAI1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ORAI1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IMD9; TAM2; ORAT1; CRACM1; TMEM142A

Catalog #: A10953

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84876

Isotype: IgG

Swiss Prot: Q96D31

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a membrane calcium channel subunit that is activated by the calcium sensor STIM1 when calcium stores are depleted. This type of channel is the primary way for calcium influx into T-cells. Defects in this gene are a cause of immune dysfunction with T-cell inactivation due to calcium entry defect type 1 (IDTICED1).

To place an order, please [Click HERE](#).