

LYZ

Reactivity:Human Mouse

Tested applications:WB IHC ICC IF IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50

Calculated MW:17 kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LYZ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LZM; LYZF1;

Catalog #:A10972

Antibody Type:

Monoclonal Antibody

Species:Rabbit

Gene ID:4069

Isotype:IgG

Swiss Prot:P61626

Purity:Affinity purification

For research use only.

Background:

This gene encodes human lysozyme, whose natural substrate is the bacterial cell wall peptidoglycan (cleaving the beta

To place an order, please [Click HERE](#).