

MCL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP FC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 37/29 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MCL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1/EAT;

Catalog #: A10993

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 4170

Isotype: IgG

Swiss Prot: Q07820

Purity: Affinity purification

For research use only.

Background:

This gene encodes an anti-apoptotic protein, which is a member of the Bcl-2 family. Alternative splicing results in multiple transcript variants. The longest gene product (isoform 1) enhances cell survival by inhibiting apoptosis while the alternatively spliced shorter gene products (isoform 2 and isoform 3) promote apoptosis and are death-inducing.

To place an order, please [Click HERE](#).