

## MKI67

**Reactivity:** Human Mouse

**Tested applications:** WB IHC ICC IF IP FC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100  
IP 1:20 - 1:50 FC 1:20 - 1:50

**Calculated MW:** 359 kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human MKI67

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KIA; MIB-; MIB-1; PPP1R105;

**Catalog #:** A11005

**Antibody Type:**

Monoclonal Antibody

**Species:** Rabbit

**Gene ID:** 4288

**Isotype:** IgG

**Swiss Prot:** P46013

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

*To place an order, please [Click HERE](#).*