

NGFR

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50

Calculated MW: 75 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NGFR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR;

Catalog #: A11050

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 4804

Isotype: IgG

Swiss Prot: P08138

Purity: Affinity purification

For research use only.

Background:

Nerve growth factor receptor contains an extracellular domain containing four 40-amino acid repeats with 6 cysteine residues at conserved positions followed by a serine/threonine-rich region, a single transmembrane domain, and a 155-amino acid cytoplasmic domain. The cysteine-rich region contains the nerve growth factor binding domain.

To place an order, please [Click HERE](#).