

## NRAS

---

**Reactivity:**Human Mouse Rat Zebrafish

**Tested applications:**WB IP FC

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A11059

**Antibody Type:**

Monoclonal Antibody

**Species:**Rabbit

**Gene ID:**4893

**Isotype:**IgG

**Swiss Prot:**P01111/P01112/P01116

**Purity:**Affinity purification

For research use only.